



GENEX
Laboratories

MR - 450
Medical Refrigerator

ExternalSize(W*D*H)	770*720*1980mm
Capacity	450L
Temperature Range	2°C-10°C
Temperature Accuracy	0.1°C
Display	Large screen LED display
Alarm	Audio and visual alarm for: High and low temperature ;System & Sensor failure
Refrigeration Type	Forced air refrigeration system
Refrigerant	R134a, CFC free
Condenser & Evaporator	High efficiency air-cooled condenser, material is cupreous.
Defrost	Auto defrost
Construction	Mono-assembly foaming (Rigid polyurethane insulation material)
External Material	Cold-rolled steel with anti-bacteria powder coating.
Internal Material	Cold-rolled steel with anti-bacteria powder coating.
Door	Three-ply transparent glass door with lock, Heating film preventing frosting Sealing joint strip with magnetic force
Shelves	5 Adjustable shelves made of PVC coated wire
Consumption	450W