

Multi Head Microscope GX – 510



Specification:

- Trinocular Head: 1PCS, IPD=48mm ~ 75mm
- Binocular Head: 4PCS, IPD=48mm ~ 75mm
- Eyepiece: Extra Wide Field Eyepiece WF10X / 20mm, 6PCS
- Infinity Plan Acromatic Objective: 4x, 10x, 40x(s), 100x(s, oil)
- Nosepiece: Inverse Quintuple Nosepiece
- Double Layer Mechanical Stage: 185x142 / 75x55mm
- Focusing: Coaxial Coarse & Fine Adjustment, Fine Division: 0.002mm
- Condenser: Abbe NA=0.90/0.25 Swing condenser
- Kohler Illumination: External Illumination, Halogen Lamp 12V/50W
- Pointer for teaching: Green LED Pointer, Brightness Adjustness

Standard accessory:

- Immersion oil
- Dust cover
- Instruction manual

Optional accessory:

- Achromatic Objective 20X
- 1X C-MOUNT, 0.5X C-MOUNT
- Can be equipped with Digital Camera Eyepiece
- Kohler Illumination: External Illumination, Halogen Lamp 12V/50W
- Pointer for teaching: Green LED Pointer, Brightness Adjustness