

THROMBO

Semi Automated Coagulation analyzer



SPECIFICATIONS:

Item modify: easy to add or edit items, easy to set the items parameters Calibration: single-point calibration; multipoint calibration; auto store and match with calibration curve Quality control: auto draw the quality figure Reagent: Open reagent system Reagent volume: 40µL Sample volume: 40µL Memory: Save more than 10000 samples results, and auto recovering the results when power is off during working. Repeatability: CV 5% Calculative item: S, %, PTR, INR, Fib can be get through PT, This leads to lower the cost of reagent Printer: Inner thermal printer **Test Items** PT APTT TT Fib RepT VT AT-III Protein-S Protein-C LA Hep