



## AT series

1. LCD display, 2 programs easy to choose
2. Automatically feeding water and excess water and air to be drained away.
3. Show the reason of failure on LCD directly
4. Safety Devices:
  - A. Over heat control: When inner temperature is overheated , power will be shut down automatically and alarm system
  - B. Over pressure control: it have alarm system when over pressure.
  - C. Safety door lock system
  - D. Low water level indicator :It have alarm system when low water

**we have new size AT 30 , it is 30 liter capacity**

Model	AT 14	AT 17	AT 23
Chamber Dimensions(mm)(Φ×D)	250×270	250×380	250×500
Chamber Capacity(liter)	14L	17L	23L
Overall Size(cm)(L×W×H)	48×48×42	48×56×42	8×67×42
Heater	1200W	1300W	1600W
Net Weight(KGS)	31	34.5	41
Gross Weight(KGS)	38	43	50
Packing Size(cm)	68×60×50	73×60×50	87×60×50
Minimum Desk Size(cm)	48×35	48×40	48×50
Power Supply	230V,50/60HZ,110V,126V		
Sterilize Temperature	121°C/134°C		
Sterilize Time	Please refer to different programs.		
Stainless steeling instrument container:	2 trays		