

Models: Labfuge12/24, Labfuge12SW
BENCHTOP LOW-SPEED CENTRIFUGE



12×20ml



24×15ml



12×15ml

Technical Parameters

| | | | |
|------------------|--------------------------------|-----------|--------|
| Model | LABFUGE | | |
| Max Speed | 4000r/min | | |
| Speed Accuracy | ±30rpm | | |
| NO.1 Swing Rotor | 12×15ml | 4000r/min | 2680×g |
| NO.2 Angle Rotor | 24×15ml | 4000r/min | 2680×g |
| NO.3 Angle Rotor | 12×20ml | 4000r/min | 2680×g |
| Motor | Brushless DC Motor | | |
| Timer Range | 0~99min | | |
| Noise | ≤50dB(A) | | |
| Power Supply | DC220V &110V 50Hz 2A | | |
| Dimension | Digital : 460×360×270mm(L×W×H) | | |
| Weight | Digital : 19.5kg | | |