



Rotary Microtome

Histo – 05

Product Description

- © It has a clipper –built case, the user can put blade, wax block, or other consumable on the top of case. It is easy.
- © It does not need to be lubricated and maintained frequently.
- © The handwheel on the right can be locked in any position, and the handle can be locked in the highest position; it is double safe for changing blade or tissue block.
- © The specimen orienting system: axes X and Y 8°, axis 180°.
- © It adopts a special blade holder; it is safer and more convenient to use disposable blades directly. It also extends the user life of disposable blades.

Technical Data

- © Slice thickness range: 0- 60um
 - Setting values : from 0-2 um in 0.5 um –increments
 - Form 2-10um in 1 um-increments
 - Form 10-20 um in 2 um-increments
 - Form 20-60um in 5 um-increments
- © Horizontal specimen stroke :28 mm
- © Vertical specimen stroke : 60mm
- © Slip precision: $\pm 5\%$
- © Maximal slice section: 50 × 45mm